



CLINNOVA Ultra 2



The CLINNOVA Ultra 2 is designed for ultrasound, sonography and scanning procedures but is extremely versatile, so could be used for multiple disciplines. The couch includes fixing points for optional Lithotomy Stirrups or Adjustable Leg Supports.

Models

See Page 26-27 for available models



Model NV9576-PRM
Colour Pepper

Standard Features

- Strong & sturdy design, supporting & lifting patients up to a maximum weight of 265Kg (583lbs)
- Gas assisted drop end foot section, adjustable from 0° to -84°
- Electric backrest section, adjustable from 0° to +90°
- Anti-trap gas assisted sections prevent risk of crushing and harm to patient, user or furniture
- Electric height elevation, raises from 46cm to 95cm in 15 seconds
- Variable height range or individual working heights
- Choice of independent, retractable or central locking castors
- Purpose-made precision bearings ensure quiet, reliable operation and increase product life expectancy
- Frame design enables use of mobile hoist
- Includes a hand switch as standard
- Laser cut fixing points to easily retrofit approved accessories
- Anti-microbial paint additive on framework prevents growth of bacteria, fungi and viruses
- 70cm wide Easy Clean Upholstery ideal for infection prevention
- Choice of Standard or Standard+ upholstery colours

Optional Accessories

- 6207** Double foot switch
- 6044** Lithotomy Stirrups
- 6045** Gynaecology Leg Supports
- 6006** Upper limb support
- 6007** Removable head/knee support
- 6057** Battery back-up facility
- 6054** Telescopic IV pole



Paper Roll Holder
Fitted as Standard



Model NV9676-CLS
Colour Sky Blue

CLINNOVA Ultra 2 - Models Available



NV9476-CLS

Electric Height Elevation
Electric Backrest
Gas Assisted Foot Section
Independently Braking Castors
Classic Base



NV9576-CLS

Electric Height Elevation
Electric Backrest
Gas Assisted Foot Section
Retractable Castors
Classic Base



NV9676-CLS

Electric Height Elevation
Electric Backrest
Gas Assisted Foot Section
Central Locking Castors
Classic Base



NV9476-PRM

Electric Height Elevation
Electric Backrest
Gas Assisted Foot Section
Independently Braking Castors
Premium Base



NV9576-PRM

Electric Height Elevation
Electric Backrest
Gas Assisted Foot Section
Retractable Castors
Premium Base



NV9676-PRM

Electric Height Elevation
Electric Backrest
Gas Assisted Foot Section
Central Locking Castors
Premium Base