

Height range 50cm to 105cm

Medicare Bariatric Drop End Couch









The Bariatric Drop End Couches feature electric independently operated foot sections for effortless positioning for various lower limb procedures. The electric backrest supports the patient in the sitting position. There is a choice of models with either single or dual foot sections, all with a safe working load and lifting capacity of 325Kg (715lbs).



SM3583

Electric height, backrest and single foot


SM3584



Standard Features

- Safe working load and lifting capacity of 325Kg (715lbs)
- Retractable wheel design with adjustable foot for stability

- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- Hand switch fitted as standard



Dual Adjustable Electric Footrests on SM3584

Optional Accessories

6005 Paper roll holder

Head support cushion

Side support loops (Pair)

Side support rails (Pair)

9002P Side support rail cushions (Pair) Triple foot switch for SM3583

Quad foot switch for SM3584



Available Upholstery Colours:

