



+75°

## Therapy Deluxe Bobath Couch

The top of the range Deluxe Bobath couch has a safe working load and lifting capacity of 325Kg (715lbs). This model features electric height and backrest operation and has large central locking castors with steering facility for ease of movement. There are a wide range of specific accessories available to further enhance this model.

## Standard Features

- Safe working load and lifting capacity of 325Kg (715lbs)
- Electric height operation
- Electric head section adjustable from 0° to
- Large central locking castors with steering control
- Excellent height range from 46cm to 100cm
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- Dual hand switch fitted as standard
- Choice of Standard or Standard+ upholstery colours



Model ST4642 Shown with optional Side Support Rails and Hand Support Loops Colour Light Grey






Models

Model ST4642

Colour Sky Blue

ST4642 ST4642W



## **Optional Accessories**

6004
6057
-2WF
-EPFS
6074
9002P
6041
6012
6013

Breathing Hole & Plug Battery Back-Up Double Foot Switch Perimeter Foot Switch Side Support Rails Side Support Rail Cushions Side Support Loops Upper and Lower Limb Support Patient Support Loops



Model ST4642 Shown with optional Side Support Rails and Patient Limb Support Colour Light Grey

